

Introducing

Clarity Choice_{Phz}

- Automatic Haze and Clarity measurements
- ASTM D8148 Haze in Fuels
- Correlates to ASTM D4176 Procedure 2 Haze rating
- Instrument Haze Rating (IHR) 1-6
- Proprietary Haze Clarity Index (HCI) 50-100 Quantification
- State of the Art Optics
- Rugged insertion probe technology
- No moving parts or filters to degrade
- Optional Ruggedized Tablet - allows on site calibration
- Store, print and export hundreds of results
- Output to LIMS, PC, Network or external storage devices
- Proudly made in the USA



Available exclusively from:



Choice Analytical, Inc.

Support You Can Trust

For product and sales information contact: sales@choiceanalytical.com or 832-271-7448
www.choiceanalytical.com